

Anti-CD106 antibody (25-310 aa) (STJ11105075)
STJ11105075

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	Vcam1
Uniprot ID	VCAM1_MOUSE
Immunogen	
Immunogen Region	25-310 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 25-310 of mouse CD106. (NP_035823.3).
Immunogen Sequence	FKIEISPEYKTIAQIGDSMA LTCSTTGCEPLFSWRTQID SPLNAKVRTEGSKSVLTMEP VSFENEHSYLCTATCGSGKL ERSIHVDIYSFPKDPEIQFS GPLEVGKPVTVKCLAPDIYP VYRLEIDLFGKQQLMNRQEF SSEEMTKSLETKSLEVTFTP VIEDIGKALVCRAKLHIDQI DSTLKERETVKELQVYISPR NTTISVHPSTRLQEGGAVTM TCSSEGLPAPEIFWGRKLD



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081