

Anti-ANXA2 antibody (1-339 aa) [S5083RM] (STJ11105083)

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5083RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	302
Gene Symbol	ANXA2
Uniprot ID	ANXA2_HUMAN
Immunogen	
Immunogen Region	1-339 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-339 of human Annexin A2. (NP_004030.1).
Immunogen Sequence	MSTVHEILCKLSLEGDHSTP PSAYGSVKAYTNFDAERDAL NIETAIKTKGVDEVITVNL TNRSNAQRQDIAFAYQRRTK KELASALKSALSGHLETIVL GLLKTPAQYDASELKASMKG LGTDEDSLIEICSRNTQEL QEINRVYKEMYKTDLEKDII SDTSGDFRKLMLVALAKGRRA EDGSVIDYELIDQDARDLYD AGVKRKGTDVPKWISIMTER SVPHLQKVFDRYKSYSPYD



Western blot analysis of various lysates, using Annexin A2 Rabbit monoclonal antibody (STJ11105083) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000859) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081