

Anti-CDKN2A/p16INK4a antibody (1-156 aa) [S5116RM] (STJ11105116)

STJ11105116

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5116RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1029
Gene Symbol	CDKN2A
Uniprot ID	CDN2A_HUMAN
Immunogen	1-156 aa
Region	
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human CDKN2A/p16INK4a. (NP_000068.1).
Immunogen Sequence	MEPAAGSSMEPSADWLATAA ARGRVEEVRLLEAGALPNA PNSYGRRPQVMMMGSAARVA ELLLLHGAEPNCADPATLTR PVHDAAREGFLDTLVVLHRA GARLDVRDAWGRLPVDLAE LGHRDVARYLRAAAGGTRGS NHARIDAAEGPSDIPD



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081