

Anti-ADIPOQ antibody (19-244 aa) [S5123RM] (STJ11105123)

STJ11105123

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5123RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9370
Gene Symbol	ADIPOQ
Uniprot ID	ADIPO_HUMAN
Immunogen	
Immunogen Region	19-244 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-244 of human Adiponectin. (NP_004788.1).
Immunogen Sequence	ETTTQGPVLLPLPKGACTG WMAGIPGHPGHNGAPGRDGR DGTGPEKGEKGDPLIGPKG DIGETGVPGAEGPRGFPGIQ GRKGPEGEGAYYRSAFSVG LETYVTIPNMPIRFTKIFYN QQNHYDGGSTGKFHCNIPGLY YFAYHITVYMKDVKVSLFKK DKAMLFTYDQYQENNVDQAS GSVLLHLEVGDQVWLQVYGE GERNGLYADNDNDSTFTGFL LYHDTN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081