

Anti-GRP antibody (STJ11105168)

STJ11105168

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2922
Gene Symbol	GRP
Uniprot ID	GRP_HUMAN
Immunogen	MRGRELPVLALLVLCAPR GRAVPLPAGGGTVLTKMYPR GNHWAVGHLMGKKSTGESSS VSERGLKQQLREYIRWEEA ARNLLGLIEAKENRNHQPPQ PKALGNQQPSWDSSESSNFK DVGSKGKVGRLSAPGSQREG RNPQLNQQ
Immunogen Region	
Specificity	Recombinant protein of human progastrin-releasing peptide.
Immunogen Sequence	MRGRELPVLALLVLCAPR GRAVPLPAGGGTVLTKMYPR GNHWAVGHLMGKKSTGESSS VSERGLKQQLREYIRWEEA ARNLLGLIEAKENRNHQPPQ PKALGNQQPSWDSSESSNFK DVGSKGKVGRLSAPGSQREG RNPQLNQQ



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081