

Anti-Hmox1 antibody (1-266) (STJ11105174)
STJ11105174

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	15368
Gene Symbol	Hmox1
Uniprot ID	HMOX1_MOUSE
Immunogen	
Immunogen Region	1-266
Specificity	Recombinant Protein corresponding to a sequence within amino acids 1-266 of mouse Hmox1 (NP_034572.1).
Immunogen Sequence	ERPQPDSPMPQDLSEALKEAT KEVHIQAENAEFMKNFQKGQ VSREGFKLVMSLYHIYTL EEEIERNKQNPVYAPLYFPE ELHRRAALEQDMAFWYGPWH QEIIPTPATQHYVKRLHEV GRTHPELLVAHAYTRYLGDL SGGQVLKKIAQKAMALPSSG EGLAFFTFPNIDSPTKFKQL YRARMNTLEMTPEVKHRVTE EAKTAFLLNIELFEELQVML TEEHKDQSPSQMASLRQRP



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081