

Anti-CD119 antibody (26-250 aa) (STJ11105176)

STJ11105176

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	15979
Gene Symbol	Ifngr1
Uniprot ID	INGR1_MOUSE
Immunogen	
Immunogen Region	26-250 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 26-250 of mouse CD119. (NP_034641.1).
Immunogen Sequence	ALTSTEDPEPPSPVPPTNVL IKSYNLNPVVCWEYQNMSQT PIFTVQVKVYSGSWTDSCNTN ISDHCCNIYEQIMYPDVSAW ARVKAKVGQKESDYARSKEF LMCLKGKVGPPGLEIRRKKE EQLSVLVFHPVEVVNGESQG TMFGDGSTCYTFDYTVYVEH NRSGEILHTKHTVEKEECNE TLCELNISVSTLDSRYCISV DGISSFQVQRTEKSKDVCIP PFHDD



Western blot analysis of EL4, using CD119 Rabbit polyclonal antibody (STJ11105176) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081