

Anti-PRKCE antibody (488-737 aa) [S5207RM] (STJ11105207)

STJ11105207

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5207RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5581
Gene Symbol	PRKCE
Uniprot ID	KPCE_HUMAN
Immunogen	
Immunogen Region	488-737 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 488-737 of human PKC Epsilon (NP_005391.1). YVNGDLMFQIORSRKDEP RSRFYAAEVTSALMFLHQHG VIYRDLKLDNILLDAEGHCK LADFGMCKEGILNGVTTTTF CGTPDYIAPEILQELEYGPS VDWWALGVLMYEMMAGQPPF EADNEDDLFESILHDDVLYP VWLSKEAVSILKAFMTKNPH KRLGCVASQNGEDAIKQHPF FKEIDWVLEQKKIKPPFKP RIKTKRDVNNFDQDFTREEP VLTLYDEAIVKQINQEEFK
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081