

Anti-AP2B1 antibody (752-951 aa) [S5219RM] (STJ11105219)

STJ11105219

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5219RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	163
Gene Symbol	AP2B1
Uniprot ID	AP2B1_HUMAN
Immunogen	
Immunogen Region	752-951 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 752-951 of human AP2B1. (NP_001025177.1).
Immunogen Sequence	MEMNFTNKALQHMTDFAIQF NKNSFGVIPSTPLAIHTPLM PNQSIDVSLPLNTLGPVMKM EPLNNLQVAVKNNIDVFYFS CLIPNLNVLVEDGKMERQVF LATWKDIPNENELQFIKEC HLNADTVSSKLQNNNVYTIA KRNVEGQDMLYQSLKLTNGI WILAEELRIQPGNPNTLSLK CRAPEVSQYIYQVYDSILKN



Western blot analysis of various lysates, using AP2B1 Rabbit monoclonal antibody (STJ11105219) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081