

Anti-HYAL1 antibody (22-435) [S5225RM] (STJ11105225)

STJ11105225

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5225RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3373
Gene Symbol	HYAL1
Uniprot ID	HYAL1_HUMAN
Immunogen	
Immunogen Region	22-435
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 22-435 of human HYAL1 (NP_149349.2).
Immunogen Sequence	FRGPLLPNRPFTTVWNANTQ WCLERHGVDDVSVFDVVAN PGQTFRGPDMTIFYSSQLGT YPYTPTGEPVFGGLPONAS LIAHLARTFQDILAAIPAPD FSGLAVIDWEAWRPRAWFNW DTKDIYRQSRALVQAQHPD WPAPQVEAVAQDQFQGAARA WMAGTLQLGRALRPRGLWGF YGFPDCYNYDFLSPNYTGQC PSGIRAQNDQLGWLWGQSRA LYPSIYMPAVLEGTKSQM



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081