

Anti-OMA1 antibody (421-520) (STJ11105230)

STJ11105230

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:5000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	115209
Gene Symbol	OMA1
Uniprot ID	OMA1_HUMAN
Immunogen	
Immunogen	421-520
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 421-520 of human OMA1 (NP_660286.1).
Immunogen	EFVDSLHGQPKMPEWLSTHP SHGNRVEYLDRLIPQALKIR EMCNCPPLSNPDPRLLFKLS TKHFLEESEKEDLNITKKQK
Sequence	MDTLPIQKQEQIPLTYIVEK



Western blot analysis of various lysates, using OMA1 Rabbit polyclonal antibody (STJ11105230) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081