

Anti-NCOA7 antibody (777-942) [S5263RM] (STJ11105263)

STJ11105263

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5263RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	135112
Gene Symbol	NCOA7
Uniprot ID	NCOA7_HUMAN
Immunogen	
Immunogen Region	777-942
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 777-942 of human NCOA7 (NP_861447.3). LRPHSALLENMHIEQLARRL PARVQGYPWRLAYSTLEHGT SLKTLRYKSASLDSPVLLVI KDMDNQIFGAYATHPFKFS
Immunogen Sequence	HHYGTGETFLYTFSPHFKVF KWSGENSYFINGDISLELG GGGGRFGLWLDADLYHGRSN SCSTFNNDILSKKEDFIVQD LEVWAFD



Western blot analysis of Mouse brain, using NCOA7 Rabbit monoclonal antibody (STJ11105263) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081