

Anti-KYNU antibody (1-465) [S5268RM] (STJ11105268)
STJ11105268

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5268RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8942
Gene Symbol	KYNU
Uniprot ID	KYNU_HUMAN
Immunogen	
Immunogen Region	1-465
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-465 of human KYNU (NP_003928.1). MEPSSLELPADTVQRIAAEL KCHPTDERVALHLDDEEDKLR HFRECFYIPKIQLPPVDLS LVNKDENAIYFLGNSLGLQP
Immunogen Sequence	KMKVKTYLEEELDKWAKIAAY GHEVGKRPWITGDESIVGLM KDIVGANEKEIALMNALTVN LHLLMLSFFKPTPKRYKILL EAKAFPSDHYAIESQLQLHG LNIEESMRMIKPREGEETLR IEDILEVIEKEGDSIAVILF SGVHFYTGQHFNPATKA



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081