

**Anti-MYH9 antibody (1860-1960) (STJ11105295)**  
STJ11105295

**GENERAL INFORMATION**


|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | IF/ICC/ELISA       |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human/Mouse        |

**PRODUCT PROPERTIES**

|                            |                                                                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                                                                  |
| <b>Clone ID</b>            |                                                                                                                                                             |
| <b>Concentration</b>       | Lot specific                                                                                                                                                |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                |
| <b>Purification</b>        | Affinity purification                                                                                                                                       |
| <b>Dilution Range</b>      | IF/ICC:1:50-1:200<br>ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.                                                                                                            |
| <b>Isotype</b>             | IgG                                                                                                                                                         |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                              |

**TARGET INFORMATION**

|                           |                                                                                                           |
|---------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 4627                                                                                                      |
| <b>Gene Symbol</b>        | MYH9                                                                                                      |
| <b>Uniprot ID</b>         | MYH9_HUMAN                                                                                                |
| <b>Immunogen</b>          |                                                                                                           |
| <b>Immunogen Region</b>   | 1860-1960                                                                                                 |
| <b>Specificity</b>        | Recombinant Protein corresponding to a sequence within amino acids 1860-1960 of human MYH9 (NP_002464.1). |
| <b>Immunogen Sequence</b> | QYKDQADKASTRLKQLKRQL EEAEEEAQRANASRRKLQRE LEDATETADAMNREVSSLKN KLRRGDLFPFVPRMARKGA GDGSDEEVDGKADGAEAKPA E |



Immunofluorescence analysis of NIH/3T3 cells using MYH9 Rabbit polyclonal antibody (STJ11105295) at a dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.

Immunofluorescence analysis of HeLa cells using MYH9 Rabbit polyclonal antibody (STJ11105295) at a dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081