

Anti-ICOSLG antibody (19-256) (STJ11105325)

STJ11105325

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB: 1:500-1:1000
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23308
Gene Symbol	ICOSLG
Uniprot ID	ICOSL_HUMAN
Immunogen	
Immunogen Region	19-256
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-256 of human ICOSL (NP_056074.1).
Immunogen Sequence	DTQKEVVRAMVGSDELSCA CPEGSRFDLNDVYVWQTSE SKTVVTHIPQNSSENVDS RYRNRALMSPAGMLRGDFSL RLFNVTPQDEQKFHCLVLSQ SLGFQEVLSVEVTLHVAANF SVPVVSAPHSPSQDELFTFC TSINGYPRPNVYWINKTDNS LLDQALQNDTVFLNMRGLYD VSVLRIARTPSVNIGCCIE NVLLQQNLTVGSQTGNDIGE RDKITENPVSTGEKNAAT



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081