

Anti-AKT1/AKT2/AKT3 antibody (1-123) (STJ11105345)
STJ11105345

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	207 10000 208 AKT1 AKT3 AKT1_HUMAN AKT3_HUMAN AKT2_HUMAN
Immunogen	
Immunogen Region	1-123
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-123 of human AKT1/AKT2/AKT3 (NP_005154.2).
Immunogen Sequence	MSDVAIVKEGWLHKRGEYIK TWRPRYFLLKNDGTFIGYKE RPQDQVDQREAPLNNFSVAQC QLMKTERPRPNTFIIRCLQW TTVIERTFHVETPEEREWEW TAIQTVDGLKKQEEEMDF RSG



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081