

Anti-KISS1 antibody (20-138) (STJ11105348)

STJ11105348

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3814
Gene Symbol	KISS1
Uniprot ID	KISS1_HUMAN
Immunogen	
Immunogen Region	20-138
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-138 of human KISS1 (NP_002247.3).
Immunogen Sequence	EPLEKVASVGNSRPTGQQL SLGLAPGEQSLPCTERKPA ATARLSRRGTSLSPPESSG SPQQPGLSAPHSRQIPAPQG AVLVQREKDLPNYNWNSFGL RFGKREAAPGNHGRSAGRG



Immunohistochemistry analysis of paraffin-embedded human placenta using KISS1 Rabbit polyclonal antibody (STJ11105348) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer, pH 6.0 before commencing with immunohistochemistry staining protocol.