

Anti-IL15 antibody (49-162) (STJ11105376)

STJ11105376

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:5000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3600
Gene Symbol	IL15
Uniprot ID	IL15_HUMAN
Immunogen	
Immunogen Region	49-162
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 49-162 of human IL15 (NP_000576.1).
Immunogen Sequence	NWVNVISDLKKIEDLIQSMH IDATLYTESDVHPCKVTAM KCFLLELQVISLESQDASIH DTVENLIILANNSLSSNGNV TESGCKECEEELEKNIKEFL QSFVHIVQMFINTS



Western blot analysis of Mouse heart, using IL15 Rabbit polyclonal antibody (STJ11105376) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081