

Anti-TNF RII/TNFRSF1B/CD120b antibody (23-258) (STJ11105379)
STJ11105379

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	21938
Gene Symbol	Tnfrsf1b
Uniprot ID	TNR1B_MOUSE
Immunogen	
Immunogen Region	23-258
Specificity	Recombinant Protein corresponding to a sequence within amino acids 23-258 of mouse TNF RII/TNFRSF1B/CD120b (NP_035740.2).
Immunogen Sequence	VPAQVVLTYPKPEPGYECQI SQEYYDRKAQMCCAKCPPGQ YVKHFCNKTSDTVCADEAS MYTQVWNQFRTCLSCSSSCT TDQVEIRACTKQQRVCACE AGRYCALKTHSGSCRQCMLR SKCGPGFGVASSRAPNGNVL CKACAPGTFSDDTTSSTDVCR PHRICSLAIPGNASTDAVC APESPTLSAIPRTLYVSQPE PTRSQPLDQEPGPSQTPSIL TSLGSTPIIEQSTKGG



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081