

Anti-DCAF1 antibody (1058-1396) (STJ11105384)

STJ11105384

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9730
Gene Symbol	DCAF1
Uniprot ID	DCAF1_HUMAN
Immunogen	
Immunogen Region	1058-1396
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1058-1396 of human VPRBP (NP_055518.1).
Immunogen Sequence	ASFPKYGGVDGGCFDRHLIF SRFRPISVFREANEDESGFT CCAFSARERFLMLGTCTGQL KLYNVFSGQEEASYNCHNSA ITHLEPSRDGSLLLTSATWS QPLSALWGMKSVFDMKHSFT EDHYVEFSKHSQDRVIGTKG DIAHIYDIQTGNKLLTLFNP DLANNYKRNCATFNPTDDL V LNDGVLWDVRSQAQIHKFDK FNMNISGVFHPNGLEVIINT EIWDLRTFHLHTVPALDQ



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081