

Anti-CCAR2 antibody (1-100) (STJ11105401)
STJ11105401

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:5000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	57805
Gene Symbol	CCAR2
Uniprot ID	CCAR2_HUMAN
Immunogen	
Immunogen Region	1-100
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CCAR2/DBC1 (NP_066997.3).
Immunogen Sequence	MSQFKRQRINPLPGGRNFSG TASTSLLGPPPGLLTPPVAT ELSQNRHLQGGEKQRVFTG IVTSLHDYFGWDEEVFFQL SVVKGRLPQLGKELVKAAY




Immunohistochemistry analysis of CCAR2/DBC1 in paraffin-embedded human colon tissue using CCAR2/DBC1 Rabbit polyclonal antibody (STJ11105401) at a dilution of 1:200 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to immunohistochemistry staining.



Immunohistochemistry analysis of CCAR2/DBC1 in paraffin-embedded human liver cancer tissue using CCAR2/DBC1 Rabbit polyclonal antibody (STJ11105401) at a dilution of 1:200 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to immunohistochemistry staining.



Immunohistochemistry analysis of CCAR2/DBC1 in paraffin-embedded rat brain tissue using CCAR2/DBC1 Rabbit polyclonal antibody (STJ11105401) at a dilution of 1:200 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to immunohistochemistry staining.



Immunohistochemistry analysis of CCAR2/DBC1 in paraffin-embedded mouse brain tissue using CCAR2/DBC1 Rabbit polyclonal antibody (STJ11105401) at a dilution of 1:200 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to immunohistochemistry staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081