

Anti-VSIR antibody (33-194) (STJ11105402)

STJ11105402

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64115
Gene Symbol	VSIR
Uniprot ID	VISTA_HUMAN
Immunogen	
Immunogen	33-194
Region	
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 33-194 of human VISTA/B7-H5/GI24 (NP_071436.1).
Immunogen Sequence	FKVATPYSLYVCPEGQNVTL TCRLGPDVKGHDVTFYKTV YRSSRGEVQTCSERRPIRNL TFQDLHLHHGGHQAANTSHD LAQRHGLASDHHGNFSIT MRNLTLLDSGLYCCLVVEIR HHHSEHRVHGAMELQVQTGK DAPSNVCVYPSSQDSENIT AA



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081