

Anti-SLC25A1 antibody (14-124) (STJ11105447)

STJ11105447

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:5000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6576
Gene Symbol	SLC25A1
Uniprot ID	TXTP_HUMAN
Immunogen	
Immunogen Region	14-124
Specificity	A synthetic peptide corresponding to a sequence within amino acids 14-124 of human SLC25A1 (NP_005975.1).
Immunogen Sequence	APASGKAKLTHPGKAILAGG LAGGIEICITFPTEYVKQL QLDERSHPPRYRGIGDCVRQ TVRSHGVLGLYRGLSLLYG SIPKAAVRFGMFEFLSNHMR DAQGRDLSTRG



Western blot analysis of various lysates using SLC25A1 Rabbit polyclonal antibody (STJ11105447) at 1:1000 dilution. Secondary antibody:HRP Goat Anti-Rabbit IgG (H+L) (STJUS000656) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection:ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081