

Anti-FABP1 antibody (1-100) (STJ11105450)

STJ11105450

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/IC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IF/IC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2168
Gene Symbol	FABP1
Uniprot ID	FABPL_HUMAN
Immunogen	
Immunogen Region	1-100
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FABP1 (NP_001434.1).
Immunogen Sequence	MSFSGKYQLSQENFEAFMK AIGLPEELIQKGKDIKGVSE IVQNGKHFKFTITAGSKVIQ NEFTVGEECELETMTGEKVK TVVQLEGNKLVTTFKNIKS



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081