

Anti-MRPL20 antibody (50-149) (STJ11105462)

STJ11105462

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:5000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55052
Gene Symbol	MRPL20
Uniprot ID	RM20_HUMAN
Immunogen	
Immunogen Region	50-149
Specificity	A synthetic peptide corresponding to a sequence within amino acids 50-149 of human MRPL20 (NP_060441.2).
Immunogen Sequence	FKVCKTARYLKKNMRTLWI NRITAASQEHGLKYPALIGN LVKCQVELNRKVLADLAIE PKTFKSLAALASRRRHEGFA AALGDGKEPEGIFSRVQYH



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081