

Anti-MRPL22 antibody (59-206) (STJ11105466)

STJ11105466

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	29093
Gene Symbol	MRPL22
Uniprot ID	RM22_HUMAN
Immunogen	
Immunogen Region	59-206
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 59-206 of human MRPL22 (NP_054899.2). GEPRRPAEIIYHCRRIKYSK DKMWYLAKLRGMSIDQALA QLEFNDKKGAKIIEVLLEA QDMAVRDHNVEFRSNLYIAE
Immunogen Sequence	STSGRGQCLKRIRYHGRGRF GIMEKVYCHYFVKLVEGPPP PPEPPKTAVAHAKEYIQQLR SRTIVHTL



Western blot analysis of lysates from 293T cells using Anti-MRPL22 Rabbit polyclonal antibody (STJ11105466) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25A Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081