

Anti-PPARA antibody (201-300) (STJ11105478)

STJ11105478

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5465
Gene Symbol	PPARA
Uniprot ID	PPARA_HUMAN
Immunogen	
Immunogen Region	201-300
Specificity	A synthetic peptide corresponding to a sequence within amino acids 201-300 of human PPARA (NP_001001928.1).
Immunogen Sequence	ADLKSLAKRIYEAYLKNFNM NKVKARVILSGKASNNPPFV IHDMETLCMAEKTIVAKLVA NGIQNKAEVRIHFCCQCTS
Sequence	VETVTELTFAKAIPGFANL



Western blot analysis of lysates from Mouse kidney using PPARA Rabbit polyclonal antibody (STJ11105478) at 1:700 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081