

Anti-PPP2R1B antibody (502-601) (STJ11105479)

STJ11105479

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5519
Gene Symbol	PPP2R1B
Uniprot ID	2AAB_HUMAN
Immunogen	
Immunogen Region	502-601
Specificity	A synthetic peptide corresponding to a sequence within amino acids 502-601 of human PP2A-A Beta/PR65 Beta/PPP2R1B (NP_002707.3).
Immunogen Sequence	ANDPNYLHRMTTLFCINALS EACGQEITTKQMLPIVLKMA GDQVANVRFNVAKSLQKIGP ILDTNALQGEVKPVLQKLQG DEDMDVKYFAQEASVLALA



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081