

Anti-PTPRC antibody (27-120) (STJ11105481)
STJ11105481

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5788
Gene Symbol	PTPRC
Uniprot ID	PTPRC_HUMAN
Immunogen	
Immunogen Region	27-120
Specificity	A synthetic peptide corresponding to a sequence within amino acids 27-120 of human CD45RO (NP_563578.2).
Immunogen Sequence	SPTSPSTGLTTAKMPVSVLS SDPLPTHITTAFFSPASTFERE NDFSETTSLSPDNTSTQVS PDSLNDNASAFNTTGVSSVQT
Sequence	PHLPTHADSQTPSA



Immunofluorescence analysis of rat spleen using Human CD45RO Rabbit polyclonal antibody (STJ11105481) at a dilution of 1:25 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining. Perform high pressure antigen retrieval with 0.01 M citrate buffer (pH 6.0) prior to IF staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081