

Anti-PTPRC antibody (10-110) (STJ11105482)

STJ11105482

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC: 1:50-1:200
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5788
Gene Symbol	PTPRC
Uniprot ID	PTPRC_HUMAN
Immunogen	
Immunogen Region	10-110
Specificity	A synthetic peptide corresponding to a sequence within amino acids 10-110 of human CD45RA (P08575-8).
Immunogen Sequence	LAFGFAFLDTEVFVTGQSPT PSPTGLTTAKMPSPVPLSSDP LPTHHTAFSPASTFERENDF SETTTSLSPDNTSTQVSPDS LDNASAFNTTDAYLNASETT T



Immunofluorescence analysis of paraffin-embedded human tonsil using Human CD45RA Rabbit polyclonal antibody (STJ11105482) at dilution of 1:50 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining. Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IF staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081