

Anti-IFI16 antibody (500-785) (STJ11105507)
STJ11105507

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/FC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	FC:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3428
Gene Symbol	IFI16
Uniprot ID	IF16_HUMAN
Immunogen	
Immunogen Region	500-785
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 500-785 of human IFI16 (NP_001351796.1).
Immunogen Sequence	KSEDTISKMNDFMRMQILKE GSHFPGPFMTSIGPAESHPH TPQMPPSTPSSSFLTLKPR LKTEPEEVSIEDSAQSDLKE VMVLNATESFVYEPKEQKKM FHATVATENEVFRVKVFNID LKEKFTPKKIIAIAANYVCRN GFLEVYPFTLVADVNADRNM EIPKGLIRSASVTPKINQLC SQTGGSFVNGVFEVHKKNVR GEFTYYEIQDNTGKMEVWH GRLTTINCEEGLKLTCTF



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081