

Anti-VEGFR3/Flt-4 antibody (25-229) (STJ11105508)

STJ11105508

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14257
Gene Symbol	Flt4
Uniprot ID	VGFR3_MOUSE
Immunogen	
Immunogen Region	25-229
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 25-229 of mouse VEGFR3/Flt-4 (NP_032055.1).
Immunogen Sequence	YSMTPTPTLNITEDSYVIDTG DLSISICRGQHPLEWTPGA QEVLTGKGKDESDTRVWHDC EGTEARPYCKVLLAQTHAN NTGSYHCYYKYIKARIEGTT AASTYVFVRDFKHPFINKPD TLLVNRKDSMWVPCLVSIPL LNITLRSQSSALHPDGGQEV WDDRRGRMRVPTQLLRDALYL QCETTWDGQNFSLNLFVVIH TGNEL



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081