

Anti-CCL8 antibody (24-99) (STJ11105512)

STJ11105512

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6355
Gene Symbol	CCL8
Uniprot ID	CCL8_HUMAN
Immunogen	
Immunogen Region	24-99
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 24-99 of human CCL8/MCP-2 (NP_005614.2).
Immunogen Sequence	QPDVSVIPITCCFNVINRKI PIQRLESYTRITNIQCPKEA VIFKTKRGKEVCADPKERWV RDSMKHLDQIFQNLKP



Western blot analysis of lysates from wild type (WT) and 293T cells transfected with CCL8/MCP-2 using CCL8/MCP-2 Rabbit polyclonal antibody (STJ11105512) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081