

Anti-GDF15 antibody (197-308) (STJ11105532)

STJ11105532

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9518
Gene Symbol	GDF15
Uniprot ID	GDF15_HUMAN
Immunogen	
Immunogen Region	197-308
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 197-308 of human GDF15 (NP_004855.2).
Immunogen Sequence	ARNGDHCPLGPGRCRLHTV RASLEDLGWADWVLSPREVQ VTCIGACPSQFRAANMHAQ IKTSLHRLKPDTVPAPCCVP ASYNPMVLIQKTDGTGVSQQT YDDLLAKDCHCI



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081