

Anti-E-Selectin/CD62E antibody (22-557) (STJ11105535)

STJ11105535

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	20339
Gene Symbol	Sele
Uniprot ID	LYAM2_MOUSE
Immunogen	
Immunogen Region	22-557
Specificity	Recombinant Protein corresponding to a sequence within amino acids 22-557 of mouse E-Selectin/CD62E (NP_035475.2).
Immunogen Sequence	WYYNASSSELMTYDEASAYCQ RDYTHLVAIQNKEEINYLNS NLKHSPSYWIGIRKVNWW IWVGTGKPLTEEAQNWAPGE PNNKQRNEDCVEIYIQRKD SGMWNDERCNKKKALCYTA SCTNASCSGHGCETINSY TCKCHPGFLGPNCQAVTCK PQEHPDYGSLNCSHPFGPFS YNSSCSFGCKRGYLPSSMET TVRCTSSGEWSAPAPACHVV ECEALTHPAHGIRKCSSNP



Western blot analysis of lysates from wild type (WT) and 293T cells transfected with E-Selectin/CD62E using A E-Selectin/CD62E Rabbit polyclonal antibody (STJ11105535) at a 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081