

Anti-DHX30 antibody (336-501) (STJ11105560)

STJ11105560

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	22907
Gene Symbol	DHX30
Uniprot ID	DHX30_HUMAN
Immunogen	
Immunogen Region	336-501
Specificity	Recombinant Protein corresponding to a sequence within amino acids 336-501 of human DHX30 (NP_619520.1).
Immunogen Sequence	RELGETQRRPCTIQVPEPIL RKIETFLNHYPVESSWIAPE LRLQSDDILPLGKDSGPLSD PITGKPYVPLLEAEVRLSQ SLLLELWRRRGPVWQEAPQLP VDPHRDTILNAIEQHPVVI SGGTGCCKTTRIPQLLLERY VTEGRGARCNVITQPRRIS AVSVAQ



Western blot analysis of various lysates using DHX30 Rabbit polyclonal antibody (STJ11105560) at 1:1000 dilution. Secondary antibody:HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection:ECL Basic Kit. Exposure time: 45s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081