

Anti-BRAF antibody (507-606) (STJ11105588)

STJ11105588

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:100-1:500 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	673
Gene Symbol	BRAF
Uniprot ID	BRAF_HUMAN
Immunogen	
Immunogen Region	507-606
Specificity	A synthetic peptide corresponding to a sequence within amino acids 507-606 of B-Raf V600E (NP_004324.2).
Immunogen Sequence	KTRHVNILLFMGYSTKPQLA IVTQWCEGSSLYHHLHIET KFEMIKLIDIA RQTAQGMDY LHAKSIIHRDLKSNINFLHE DLTVKIGDFGLATEKSRWSG