

Anti-CD71 antibody (89-763) (STJ11105658)

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	22042
Gene Symbol	Tfrc
Uniprot ID	TFR1_MOUSE
Immunogen	
Immunogen Region	89-763
Specificity	Recombinant Protein corresponding to a sequence within amino acids 89-763 of mouse CD71 (NP_035768.1).
Immunogen Sequence	CKRVEQKEECVKLAETEED KSETMETEDVPTSSRLYWAD LKTLSEKLSNIEFADTIKQ LSQNTYTPREAGSQKDESLA YYIENQFHEFKSVWRDEH YVKIQVKSSIGQNMVTIVQS NGNLDPVESPEGYVAFSKPT EVSGKLVHANFGTKKDFEEL SYSVNGSLVIVRAGEITFAE KVANAQSFNAIGVLIYMDKN KFPVVEADLALFGHAHLGTG DPYTPGFPSFNHTQFPSPQ



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081