

Anti-NECTIN3 antibody (127-366) (STJ11105695)

STJ11105695

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:100-1:500
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25945
Gene Symbol	NECTIN3
Uniprot ID	NECT3_HUMAN
Immunogen	
Immunogen Region	127-366
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 127-366 of human PVRL3 (NP_001230215.1). SLNDATITLHNIGFSDSGKY ICKAVTFPLGNAQSSTTVTV LVEPTVSLIKGPDSLIDGGN ETVAACIAATGKPVAHIDW EGDLGEMESTTTSFPNETAT IISQYKLFPTFRFARGRRITC VVKHPALEKDIRYSFILDQI YAPEVSVTGYDGNWVFGKRG VNLKCNADANPPPFKSVWSR LDGQWPDGLLASDNTLHFVH PLTFNYSGVYICKVTNSLGQ RSDQKVIYISAYNSVASLN
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081