

Anti-IKBIP antibody (115-350) (STJ11105730)

STJ11105730

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	IF/CC:1:100-1:500
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	121457
Gene Symbol	IKBIP
Uniprot ID	IKIP_HUMAN
Immunogen	
Immunogen Region	115-350
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 115-350 of human IKBIP (NP_963906.1).
Immunogen Sequence	IIAHLKRLQEEINEVKTWSN RITEKQDILNNSLTLSQDI TKVDQSTTSMKDVGLKITS VKTDIRISGLVTDVISLTD SVQELENKIEKVEKNTVKNI GDLLSSSIDRTATLRKTASE NSQRINSVKKLTTELKSDFD KHTDRFLSLEGDRKVLKTV TFDNLKPKVYNLKKDFSRL EPLVNDLTLRIGRLVTDLLQ REKEIAFLSEKISNLTIVQA EIKDIKDEIAHISDMN



Immunofluorescence analysis of NIH/3T3 cells using IKBIP Rabbit polyclonal antibody (STJ11105730) at dilution of 1:500 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using IKBIP Rabbit polyclonal antibody (STJ11105730) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081