

Anti-FTSJ1 antibody (170-329) (STJ11105734)
STJ11105734

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	24140
Gene Symbol	FTSJ1
Uniprot ID	TRM7_HUMAN
Immunogen	
Immunogen Region	170-329
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 170-329 of human FTSJ1 (NP_036412.1).
Immunogen Sequence	LQVFFSSVLCAPRSSRNSS IEAFVQCQGYDPPEGFIPDL SKPLLDHSYDPDFNQLDGPT RIIVPFVTCGDLSSYDSDRS YPLDEGGSEYKYPPTQPP ISPPYQEAETLKRKGQLAKE IRPQDCPISRVDTFPQPLAA PQCHTLAPEMEDNEMSCSP



Western blot analysis of various lysates, using FTSJ1 Rabbit pAb (STJ11105734) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJ5000656) at 1:10000 dilution.
Lysates/proteins: 25 Mu g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit
Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081