

Anti-Phospho-PI3KP85alpha/P55gamma/P85Beta-Y467/Y464/Y199 antibody (STJ11105781)

STJ11105781

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat


PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5295
	5296
	8503
	PIK3R1
	PIK3R2
	P85A_HUMAN
	P85B_HUMAN
	P55G_HUMAN

Immunogen	
Immunogen Region	
Specificity	A phospho specific peptide corresponding to residues surrounding Y467 of human PIK3R1/2/3
Immunogen Sequence	



Western blot analysis of lysates from 293 cells, using Phospho-PI3KP85 Alpha/P55 Gamma/P85 Beta-Y467/Y464/Y199 Rabbit pAb (STJ11105781). Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% BSA.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081