

Anti-Phospho-NFE2L2-S40 antibody (STJ11105795)
STJ11105795

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4780
Gene Symbol	NFE2L2
Uniprot ID	NF2L2_HUMAN
Immunogen	DFSQRRKEYELEKQKKLEKE RQEQLQKEQEKAFFAQLQLD EET
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around S40 of human NRF2 (NP_006155.2).
Immunogen Sequence	DFSQRRKEYELEKQKKLEKE RQEQLQKEQEKAFFAQLQLD EET



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081