

Anti-Phospho-p38 MAPK-T180/Y182 antibody [S5844RM] (STJ11105844)

STJ11105844

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5844RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:2000-1:20000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6300 5603 5600 MAPK12 MAPK13 MK12_HUMAN MK13_HUMAN MK11_HUMAN
Immunogen	MTGYV
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around T180 & Y182 of human Phospho-p38 MAPK-T180/Y182MAPK11 (NP_002742.3).
Immunogen Sequence	MTGYV



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081