

Anti-NRF2 antibody (243-400 aa) [SRM] (STJ11105915)

STJ11105915

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/ICC:1:200-1:2000
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	4780
Gene Symbol	NFE2L2
Uniprot ID	NF2L2_HUMAN
Immunogen	
Immunogen Region	243-400 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 243-400 of human NRF2 (NP_006155.2).
Immunogen Sequence	SPHFLNAFEDSFSSILSTED PNQLTVNSLNSDATVNTDFG DEFYSAFIAEPSISNSMPSP ATLSHSLSELLNGPIDVSDL SLCKAFNQNHPESTAEFNDS DSGISLNTSPSVASPEHSVE SSSYGDTLLGLSDSEVEELD SAPGSVKQNGPKTPVHSS