

Anti-VE-Cadherin antibody (259-358 aa) [SRM] (STJ11105916)

STJ11105916

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:2000-1:4000 IHC-P:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	1003 12562
Gene Symbol	CDH5 Cdh5
Uniprot ID	CADH5_HUMAN CADH5_MOUSE
Immunogen	
Immunogen Region	259-358 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 259-358 of human VE Cadherin (NP_001786.2).
Immunogen Sequence	TQTKYTFVVPEDTRVGTSVG SLFVEDPDEPQNRMTKYSIL RGDYQDAFTIETNPAHNEGI IKPMKPLDYEQYSFIVE ATDPTIDLRYMSPAGNRAQ