

Anti-BAI1 antibody (60-350 aa) (STJ11105925)

STJ11105925

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	575
Gene Symbol	ADGRB1
Uniprot ID	AGRB1_HUMAN
Immunogen	
Immunogen Region	60-350 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 60-350 of human BAI1 (NP_001693.2).
Immunogen Sequence	VFPANASRCSWTLRNPDPRR YTYMKVAKAPVPCSGPGRV RTYQFDSFLESTRTYLGVES FDEVLRRLCDPSAPLAFQAS KQFLQMRQPPQHDGLRPR AGPPGPTDDFSVEYLVVGNR NPSRAACQMLCRWLDACLAG SRSSHPCGIMQTPCACLGGE AGGPAAGPLAPRGDVCLRDA VAGGPENCLTSLTQDRGGHG ATGGWKLWLSWGECTRDCGG GLQTRTRTCLPAPGVEGGG