

Anti-SECTM1 antibody (29-150 aa) (STJ11105927)

STJ11105927

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6398
Gene Symbol	SECTM1
Uniprot ID	SCTM1_HUMAN
Immunogen	
Immunogen Region	29-150 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 29-150 of human SECTM1 (NP_002995.1).
Immunogen Sequence	QNEGWDSPICTEGVSVSWG ENTVMSCNISNAFSHVNIKL RAHGQESAIFNEVAPGYFSR DGWQLQVQGGVAQLVIKGAR DSHAGLYMWHLVGHQRNNRQ VTLEVSGAEPQSAPDTGFWP VP