

Anti-CLPTM1L antibody (32-280 aa) (STJ11105929)

STJ11105929

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	81037
Gene Symbol	CLPTM1L
Uniprot ID	CLP1L_HUMAN
Immunogen	
Immunogen Region	32-280 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 32-280 of human CLPTM1L (NP_110409.2).
Immunogen Sequence	YTRPCSGDANCIQPYLARRP KLQLSVYTTTRSHLGAENNI DLVLNVEDFDVESKFERTVN VSVPKKTRNNGTLYAYIFLH HAGVLPWHDGKQVHLVSPLT TYMVPKPEEINLLTGESDTQ QIEAEKKPTSALDEPVSHWR PRLALNVMADNFVFDGSSLP ADVHRYMKMIQLGKTVHYLP ILFIDQLSNRVKDLMLVINRS TTELPLTVSYDKVSLGRLRF WIHMQDAVYSLQQGFSEK