

Anti-GRM7 antibody (220-410 aa) (STJ11105933)

STJ11105933

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2917
Gene Symbol	GRM7
Uniprot ID	GRM7_HUMAN
Immunogen	
Immunogen Region	220-410 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 220-410 of human GRM7 (NP_870989.1).
Immunogen Sequence	YVSTLASEGSYGEGVESFT QISKEAGGLCIAQSVRIPOE RKDRTIDFDRIKQLLDTPN SRAVVIFANDEDIKQILAAA KRADQVGHFLWVGSDSWGSK INPLHQHEDIAEGAITIQPK RATVEGFDAYFTSRTLENNR RNWVFAEYWEENFNCKLTIS GSKKEDTRKCTGQERIGKD SNYEQEGKVQF